

Cy3.5 NHS ester TEA

Chemical Properties

CAS No. :

Formula: C49H60N4O16S4

Molecular Weight:

Keep away from direct sunlight

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cy3.5 NHS ester TEA is a reactive dye employed for labeling amino groups in peptides, proteins, and oligonucleotides.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481