

## DiSulfo-Cy5 alkyne TEA

### Chemical Properties

CAS No. :

Formula: C41H56N4O7S2

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	DiSulfo-Cy5 alkyne TEA is a fluorescent dye, utilized as a reagent in click chemistry.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481