

Cy5.5 hydrazide

Chemical Properties

CAS No. :

Formula: C₄₀H₄₅ClN₄O

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Cy5.5 hydrazide is a fluorescent dye containing a hydrazide group. This compound is employed for labeling aldehydes and ketones. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481