

Recombinant endoglycoceramidase I assisted by activator II

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Recombinant endoglycoceramidase I assisted by activator II (rEGCase I assisted by Activator II) is a mixture of EGCase I and Activator II. EGCase I is a specific hydrolase ideal for comprehensive analysis of glycolipid compounds (GSLs). Activator II enhances the activity of EGCase without requiring any detergents.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481