

ST6 Sialyltransferase 6

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ST6Sialyltransferase6 (ST6GALNAC6) is a member of the sialyltransferase family, capable of modifying proteins and ceramides on the cell surface, thereby altering interactions between cells or between the cells and the extracellular matrix. ST6Sialyltransferase6 is useful in studies related to inflammatory bowel disease (IBD).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481