

Recombinant endoglycoceramidase I

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Recombinant endoglycoceramidase I (rEGCase I) is a glycosidase that facilitates the hydrolysis of β -glycosidic bonds between oligosaccharides and ceramides. Additionally, it catalyzes transglycosylation reactions, transferring the sugar moiety of glycosphingolipids to the primary hydroxyl groups of various 1-alkanols.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481