

Cyanine 5 Bisfunctional MTSEA Dye potassium

Chemical Properties

CAS No. :

Formula: C43H59KN4O12S6

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Cyanine 5 Bisfunctional MTSEA Dye potassium, part of the cyanine dye family, serves as a common fluorescent marker for biomolecules, interacting effectively with them. This dye may also intercalate with double-helical DNA, exhibiting enhanced fluorescence upon binding.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481