

ST8 alpha-2,8-Sialyltransferase 4

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	ST8 alpha-2,8-Sialyltransferase 4 (ST8SIA4) is a sialyltransferase enzyme involved in the sialylation process. It is also associated with breast cancer progression, facilitating cancer cell proliferation and invasion in a sialyltransferase-dependent manner.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481