

PP18

Chemical Properties

CAS No. :

Formula: C269H426N74O88S2

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | PP18 is an innovative OXM-based dual agonist targeting GLP-1 and glucagon receptors, exhibiting an EC50 value of 0.487 nM for hGLP-1R and 0.188 nM for hGcgR. |
| Targets(IC50) | Glucagon Receptor |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481