

## pFYN peptide

## Chemical Properties

CAS No. :

Formula: C115H183N32O37PS

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	pFYN peptide is a biotin-labeled phosphorylated peptide substrate derived from FYN (a member of the Src family kinases).
Targets(IC50)	Src

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481