

## Human Serum Albumin-AF555

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**  
**Store at -20°C**  
Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Human Serum Albumin-AF555 is a conjugate of HAS with AF555 (Ex=752 nM, Em=780 nm), suitable for studying cell interactions.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481