

K90-114TAT

Chemical Properties

CAS No. :

Formula: C205H357N77O55

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	K90-114TAT is an inhibitor of the EAG2-Kv β 2 interaction, exhibiting antitumor activity against glioblastoma.
Targets(IC50)	Potassium Channel

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481