

DBCO-PEG-SH

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight
Store at -20°C
Actual storage temperature shall be subject to the COA.

Biological Description

Description	DBCO-PEG-SH is a versatile dye used extensively in biological experiments. Dyes are critical tools that assist researchers in observing and analyzing cellular structures, tracking biomolecules, assessing cell functions, distinguishing cell types, detecting biomolecules, studying histopathology, and monitoring microorganisms. Their applications span from fundamental scientific research to clinical diagnostics. Dyes also have broad uses in traditional fields like textile dyeing and emerging areas such as functional textile processing, food coloring, and dye-sensitized solar cells.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481