

Trypsin-Chymotrypsin 1:1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Trypsin-Chymotrypsin 1:1 is a biocatalyst and a pivotal enzyme in new biocatalytic technology. Enzyme engineering focuses on improving enzyme reaction kinetics, substrate selectivity, and activity under extreme conditions (e.g., low or high pH). By introducing stimulus responsiveness into these enzyme modifications, dynamic control of activity is also achievable.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481