

Aspartate Aminotransferase, Human

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Aspartate Aminotransferase, Human is a biocatalyst and serves as a pivotal enzyme in new biotechnology. Enzyme engineering aims to improve the reaction kinetics, substrate specificity, and activity under extreme conditions (such as low or high pH). By introducing stimuli-responsive modifications to these enzymes, the dynamic regulation of activity can also be achieved.
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481