

Ac-VAD-AFC

## Chemical Properties

CAS No. :

Formula: C<sub>24</sub>H<sub>27</sub>F<sub>3</sub>N<sub>4</sub>O<sub>8</sub>

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Ac-VAD-AFC serves as a fluorescent substrate for caspases and is utilized to assess caspase activity.
Targets(IC50)	Caspase

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481