

Chimax-O

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Chimax-O (Endo-chitosanase) is an enzyme of the GH46 family that catalyzes chitosan substrates to produce chitobiose and D-glucosamine (GlcN). It acts through a random endo mechanism to generate chitosan derivatives with varying degrees of polymerization (DP) ≤ 2 .
Targets(IC50)	Endogenous Metabolite

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481