

S1H

Chemical Properties

CAS No. :

Formula: C94H141N23O27

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	S1H is an analog of human growth hormone (hGH) and acts as an antagonist to the human growth hormone receptor (hGHR). It inhibits the interaction between hGH and hGHR and prevents the phosphorylation of STAT5 in cells treated with hGH.
Targets(IC50)	GHSR,STAT

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481