


PAR4 antagonist 3

Chemical Properties

CAS No. :	3057206-40-0
Formula:	C ₂₂ H ₁₆ FN ₃ O ₅ S
Molecular Weight:	453.44
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	PAR4 antagonist3 (Compound 36) is a selective antagonist of protease-activated receptor 4 (PAR4). It exhibits antiplatelet activity with an IC ₅₀ of 26.1 nM and enhances metabolic stability in human liver microsomes, with a half-life (T _{1/2}) of 97.6 minutes.
Targets(IC ₅₀)	Protease-activated Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.2054 mL	11.0268 mL	22.0536 mL
5 mM	0.4411 mL	2.2054 mL	4.4107 mL
10 mM	0.2205 mL	1.1027 mL	2.2054 mL
50 mM	0.0441 mL	0.2205 mL	0.4411 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481