

Fluorescein-PEG4-NHS ester

Chemical Properties

CAS No. :

Formula: C₃₆H₃₇N₃O₁₃S

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Fluorescein-PEG4-NHS ester, a chromogenic probe, is synthesized through the reaction of adipic acid dihydrazide and fluorescein isothiocyanate.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481