

CXCR2 Probe 1

Chemical Properties

CAS No. :

Formula: C₂₀H₂₄FN₃O₄

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	CXCR2Probe 1 (Compound[¹⁸ F]16b) is a selective ligand for CXCR2 and serves as a radioactive tracer for PET imaging of neutrophils in inflammatory diseases.
Targets(IC50)	CXCR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481