

## Cyanine 3 Bisfunctional MTSEA Dye potassium

### Chemical Properties

CAS No. :

Formula: C41H57KN4O12S6

Molecular Weight:

Keep away from direct sunlight

Storage:

Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Cyanine 3 Bisfunctional MTSEA Dye potassium, part of the cyanine dye family, is a widely used fluorescent marker for biomolecules that interacts with them effectively. Cyanine dyes can also intercalate with double-helical DNA, exhibiting enhanced fluorescence upon binding.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481