

Visepegenatide

Chemical Properties

CAS No. : 2641081-54-9
Formula: (C₂H₄O)_nC₁₉H₂₉N₅O₆S₂
Molecular Weight:
Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	PB-119 is a PEGylated exenatide that lowers blood glucose levels by enhancing β cell function and reducing insulin resistance. It is applicable in studies related to type 2 diabetes.
Targets(IC50)	IGF-1R

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481