

## Biotin-YVAD-FMK

## Chemical Properties

CAS No. :

Formula: C32H45FN6O9S

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Biotin-YVAD-FMK is a vacuolar processing enzyme (VPE) inhibitor that suppresses 64% of VPE activity at a concentration of 5 $\mu$ M. It also acts as an irreversible substrate for caspase 1.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481