

Boc-NHCH₂-Ph-pyrimidine-NH₂

Chemical Properties

CAS No. : 3040121-15-8


Formula: C₁₆H₂₀N₄O₂

Molecular Weight: 300.36

Keep away from direct sunlight

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Boc-NHCH ₂ -Ph-pyrimidine-NH ₂ serves as the linker for PROTACXYD198.
Targets(IC ₅₀)	PROTAC Linker

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.3293 mL	16.6467 mL	33.2934 mL
5 mM	0.6659 mL	3.3293 mL	6.6587 mL
10 mM	0.3329 mL	1.6647 mL	3.3293 mL
50 mM	0.0666 mL	0.3329 mL	0.6659 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481