

FA2G2 glycan (G2F), APTS labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FA2G2 glycan (G2F), APTS labelled (FA2G2 N-linked oligosaccharide, APTS labelled; F(6) A2G2, APTS labelled) is a polymeric N-glycoprotein that serves as a versatile fluorescent linker. The resultant conjugate, simulating the antennae elements of N-glycans, offers high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481