

FA2[3]BG1 & FA2[6]BG1 glycan (G1B), 2-AA labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	The compound FA2[3]BG1 & FA2[6]BG1 glycan (G1B), 2-AA labelled is a poly N-glycoprotein functioning as a versatile fluorescent linker. This conjugate is designed with high sensitivity and specificity by mimicking the antennae elements of N-linked oligosaccharides.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481