

FA2G2 glycan (G2F), 2-AB labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	FA2G2 glycan (G2F), labeled with 2-AB (FA2G2 N-linked oligosaccharide, 2-AB labelled; F(6)A2G2, 2-AB labelled), is a polymeric N-glycan protein that serves as a versatile fluorescent linker. The resulting conjugate, through mimicking the antenna elements of N-glycans, exhibits high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481