

## M5 glycan, 2-AA labelled

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	M5 glycan (Man5), 2-AA labelled (Mannose-5 N-linked oligosaccharide, 2-AA labelled; Oligomannose 5 glycan, 2-AA labelled) is a type of N-linked glycoprotein that acts as a versatile fluorescent linker. The resulting conjugates exhibit high sensitivity and specificity by mimicking the antenna structures of N-glycans.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481