

## A2[3]G1 & A2[6]G1 glycan (G1)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	A2[3]G1 & A2[6]G1 glycan (G1) is a type of poly N-glycoprotein that serves as a multifunctional fluorescent linker. The resulting conjugate simulates the antennae elements of N-glycans, providing high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481