

A3G3S3 glycan, 2-AA labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	A3G3S3 glycan, 2-AA labelled (A3G3S3 N-linked oligosaccharide, 2-AA labelled), is a type of N-glycoprotein polymer acting as a versatile fluorescent linker (Linker). The resultant conjugates exhibit high sensitivity and specificity by mimicking the antennae components of N-glycans.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481