

Alpha-Gal standard, 2-AB labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Alpha-Gal standard, 2-AB labelled (Galactose-alpha-1,3-galactose, 2-AB labelled) is a polyglycosylated N-glycoprotein that serves as a multifunctional fluorescent linker. The resulting conjugates, which simulate the antenna components of N-glycans, exhibit high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481