

## Anti-ANGPT2 Antibody

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-ANGPT2 Antibody is a humanized monoclonal antibody expressed in CHO cells targeting Angiopoietin-2 (Ang2/ANGPT2). This product, featuring a hulgG1 heavy chain and a hukappa light chain, has a predicted molecular weight of 150 kDa. It binds specifically to Ang2 with high affinity, blocking its interaction with the Tie2 receptor. This blockade removes Ang2-mediated vascular destabilization signals, inhibiting vascular endothelial cell activation, angiogenesis, and vascular leakage.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481