

## Anti-Meflin Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Meflin Antibody is a specific monoclonal antibody targeting Meflin (also known as ISLR), a marker for mesenchymal stem cells and a subset of cancer-associated fibroblasts (CAFs). This antibody is primarily used to study the role of stromal components in tumor progression, organ fibrosis, and stem cell niche regulation, particularly in the research of stroma remodeling in scirrhous tumors like pancreatic cancer.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481