

BMS-986178

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	BMS-986178 is a fully human OX40 (CD134) agonistic monoclonal antibody. By binding to the OX40 receptor on activated T cells, it provides essential co-stimulatory signals that enhance T cell expansion, survival, and cytokine secretion, thereby achieving anti-tumor effects in various solid tumor studies by recruiting effector immune cells.
Targets(IC50)	TNF

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481