

Anti-oxLDL Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-oxLDL Antibody is a specific monoclonal antibody targeting oxidized low-density lipoprotein (oxLDL). As oxLDL is a key pro-inflammatory factor promoting atherosclerotic plaque formation and vessel wall inflammation, this antibody is frequently used to study the pathological mechanisms of cardiovascular disease, plaque stability, and the development of diagnostic assays by neutralizing oxLDL or blocking its interaction with macrophage receptors (e.g., CD36).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481