

IMA-026

## Chemical Properties

CAS No. : 1038541-65-9

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	IMA-026 represents a highly specialized humanized monoclonal antibody specifically engineered to target the human interleukin-13 cytokine with exceptional binding affinity and selectivity, effectively facilitating the measurable blockade of ligand interaction with both the IL-13R alpha 1 and IL-13R alpha 2 subunits across various preclinical experimental models to evaluate the modulation of Th2-driven inflammatory responses during strictly monitored laboratory observation periods.
Targets(IC50)	IL Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481