

Anti-PCLA Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-PCLA Antibody is a specialized monoclonal antibody target the human C-type lectin domain family 14 member A protein, inhibition of tumor-associated angiogenesis and endothelial cell migration across various preclinical experimental oncology models to suppress the progression of vascularized solid tumors.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481