

## Anti-Amyloid Beta Antibody (DLX212)

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-Amyloid Beta Antibody (DLX212) represents a highly specialized humanized monoclonal antibody specifically engineered to target various forms of Amyloid Beta (A $\beta$ ) with exceptional binding affinity. It effectively facilitates the measurable clearance of amyloid plaques across various preclinical experimental models to evaluate the reduction of neurotoxicity and the restoration of synaptic plasticity during strictly monitored laboratory observation periods and functional cognitive assays.
Targets(IC50)	Beta Amyloid

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481