

## Anti-ICAM3/CD50 Antibody

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-ICAM3/CD50 Antibody is a highly specialized monoclonal antibody target the human intercellular adhesion molecule 3, which is also known as CD50, inhibition of ligand-receptor interactions between leukocytes and endothelial cells across various preclinical experimental models to modulate immune cell adhesion and inflammatory responses.
Targets(IC50)	Integrin
In vitro	In T lymphocyte and HUVEC co-culture assays, Anti-ICAM3/CD50 Antibody inhibited leukocyte adhesion and blocked tyrosine phosphorylation of intracellular proteins [1].

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481