

## Anti-TNFRSF5/CD40 Antibody

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Anti-TNFRSF5/CD40 Antibody is a highly specialized monoclonal antibody target the human tumor necrosis factor receptor superfamily member 5, which is also widely known as CD40, modulation of co-stimulatory signaling pathways and immune cell activation across various preclinical experimental models to regulate humoral immunity and dendritic cell maturation.
Targets(IC50)	TNF
In vitro	In human B lymphocyte activation assays, Anti-TNFRSF5/CD40 Antibody (mAb G28-5) combined with IL-4 induced sustained cell proliferation and cluster formation [1].

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481