

Anti-Complement C5aR1 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-Complement C5aR1 Antibody is an antagonistic monoclonal antibody targeting C5a receptor 1 (C5aR1). It inhibits downstream chemotaxis and inflammatory responses by neutralizing C5a-receptor interaction.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481