

Anti-CD7 Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Anti-CD7 Antibody is a monoclonal antibody specifically recognizing the CD7 surface antigen, which is highly expressed in T-cell acute lymphoblastic leukemia (T-ALL) and some acute myeloid leukemias. This antibody is primarily used for developing targeted therapeutic strategies, such as antibody-drug conjugates (ADCs) or nanobody drugs, aimed at precisely clearing malignant T-lymphocytes while minimizing toxicity to normal hematopoietic stem cells.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481