

IMC-C103C

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	IMC-C103C is a highly specialized bispecific T-cell receptor fusion protein designed to target the MAGE-A4 peptide presented by HLA-A*02:01 while simultaneously engaging the CD3 complex on T lymphocytes to facilitate the measurable redirection of polyclonal effector T cells toward antigen-positive tumor cells across various preclinical experimental models.
Targets(IC50)	Integrin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481