

AMG-811

Chemical Properties

CAS No. : 1439393-01-7

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | AMG-811 is a fully humanized monoclonal antibody expressed by CHO cells, which specifically targets interferon- γ (IFN- γ) and inhibits the downstream inflammatory signaling pathway by neutralizing its binding to the receptor. AMG-811 is often used to study the immune response mediated by IFN- γ and its role in autoimmune diseases (such as systemic lupus erythematosus). |
| Targets(IC50) | COX |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481