

## Anti-RSPO1 Antibody

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature  
Store at -20°C

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	Anti-RSPO1 Antibody is a specific monoclonal antibody targeting R-spondin-1 (RSPO1), a potent activator of the Wnt/beta-catenin signaling pathway. RSPO1 enhances signaling by binding to LGR4/5/6 receptors. This antibody specifically neutralizes RSPO1 and blocks its-mediated signal activation, serving as a vital tool for studying stem cell proliferation, intestinal crypt maintenance, and the progression of Wnt-dependent tumors.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481