

Anti-RHD/CD240d Antibody

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Store at low temperature
Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-RHD/CD240d Antibody is a highly specialized monoclonal antibody target the human Rh blood group D antigen (RHD/CD240d), effectively facilitating the targeted clearance of RhD-positive erythrocytes through the activation of antibody-dependent cellular cytotoxicity and phagocytosis across various preclinical immunological models.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481