

M7 glycan, 2-AB labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	M7 glycan (Man7), 2-AB labelled (Mannose-7 N-linked oligosaccharide, 2-AB labelled; Oligomannose 7 glycan, 2-AB labelled) is a polymeric N-linked glycoprotein that serves as a versatile fluorescent linker. The resulting conjugate, by mimicking the antennae elements of N-glycans, exhibits high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481