

FA2[3]BG1 & FA2[6]BG1 glycan (G1B), 2-AB labelled

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	FA2[3]BG1 & FA2[6]BG1 glycan (G1B), 2-AB labelled, is a multifunctional fluorescent linker derived from glycan N-linked oligosaccharides, specifically G1F with bisecting GlcNAc. The resulting conjugates imitate the antennary structures of N-glycans, offering high sensitivity and specificity.
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481