

Sialylated Core 1 O-glycan (C1S(3)1), 2AB labelled

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Sialylated Core 1 O-glycan (C1S(3)1), 2AB labeled, is a multifunctional fluorescent linker comprised of a polymeric N-glycoprotein. The resulting conjugates, which emulate the antenna components of N-glycans, exhibit high sensitivity and specificity.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481